

# Next Generation SedLine<sup>®</sup> Brain Function Monitoring

More Complete Data, Featuring an Enhanced Patient State Index (PSi)



## Next Generation SedLine helps clinicians monitor the state of the brain under anaesthesia with:

- > An enhanced signal processing engine, which improves performance of the Patient State Index (PSi)
- > Enhanced PSi with less susceptibility to EMG interference and improved performance in low power EEG cases
- > Four simultaneous channels of frontal electroencephalogram (EEG) waveforms, enabling bilateral data acquisition and processing of EEG signals
- > A Multitaper Density Spectral Array (DSA), which may enhance visibility of EEG features




# Improved Patient State Index (PSi)

Next Generation SedLine features an enhanced signal processing engine which provides an enhanced Patient State Index (PSi), a processed EEG parameter related to the effect of anaesthetic agents.

## Expert Scoring of Next Generation SedLine<sup>1</sup>

EEG experts scored the improvement in PSi performance between the original SedLine PSi and Next Generation SedLine PSi.

Experts found a **25% average improvement** in Next Generation PSi performance.



*To evaluate the performance of Original PSi and Next Generation PSi, independent EEG experts reviewed validation cases with both Original PSi and Next Generation PSi (blinded to the version), along with additional clinical information (MOAAS scores, EEG waveforms, drug doses, and vital signs). Compared to the expert-assessed anaesthetic depth, an error was defined as a case when expert assessment of PSi was 'Low' or 'High' and success was defined as a case when the expert assessment of PSi was 'Good'.*




**Next Generation SedLine Brain Function Monitoring** helps clinicians monitor the state of the brain under anaesthesia with bilateral data acquisition and processing of four leads of electroencephalogram (EEG) signals.

# Reducing Electromyography (EMG) Susceptibility


EMG is a common confounding factor that can interfere with EEG signals used in brain function monitoring.<sup>2</sup> Researchers have found that EMG interference existed in up to 38% of monitored patients.<sup>3</sup>

Next Generation SedLine utilizes Masimo's Parallel Signal Processing Engines to extract a clearer EEG signal and compute a processed EEG parameter (PSi) less influenced by EMG.




This image captures a moment when Next Generation SedLine detects EMG in the two engines depicted.

The case below demonstrates Next Generation SedLine's improvement to PSi in the presence of EMG interference.<sup>1</sup>




# Improved PSi Performance in Low Power EEG

Power across all frequency bands decreases with age.<sup>4</sup> Low power EEG can present a challenge for conventional brain function monitors.<sup>2</sup>




The subjects above were administered Propofol and were in a comparable anaesthetic state.<sup>5</sup>

When computing PSi, Next Generation SedLine uses adaptive signal processing with band-independent features to search for EEG features across many frequency bands and offer improved PSi performance in cases of lower power EEG.




The case below demonstrates Next Generation SedLine's improvement to PSi in low power EEG.<sup>1</sup>



# Multitaper Density Spectral Array (DSA)


Next Generation SedLine offers clinicians the flexibility of choosing to display either an enhanced Multitaper Density Spectral Array (DSA) or a standard Hanning DSA. The DSA contains left and right spectrograms representing the power of the EEG on both sides of the brain.

When using a Multitaper DSA, EEG data are transformed into the frequency domain, which may provide a better display of EEG features.



# Next Generation SedLine on Root®

The Next Generation SedLine module easily plugs into the Root patient monitoring platform via Masimo Open Connect® (MOC-9®) ports. Root's customisable, easily-interpretable display offers multiple views of brain monitoring information expanding visibility in the operating room and intensive care unit.



## A More Complete Picture of the Brain

Next Generation SedLine can be used simultaneously with O3® Regional Oximetry on the Root platform for a more complete picture of the brain.




## RD SedLine™ EEG Sensor

- > Soft foam pads for comfortable application on a patient's forehead
- > Adhesive anchors ensure secure sensor placement for optimal signal quality
- > Pre-filled gel electrodes help streamline sensor application workflows
- > Durable, heat-laminated flex cables allow flexible placement and repositioning



Four active EEG leads on the RD SedLine sensor collect data in the frontal lobe



- < Application graphics for O3 regional oximetry sensor placement simplify simultaneous application of both monitoring technologies

## SedLine Specifications

PHYSICAL CHARACTERISTICS	ENVIRONMENTAL
<b>Module Physical Dimensions</b> Width ..... 1.3 in (3.3 cm) Length ..... 4.0 in (10.2 cm) Thickness ..... 0.8 in (2.0 cm)	<b>Module Operating Conditions</b> Operating Temperature ..... 41-104°F (5-40°C) Operational Humidity ..... 15-95%, non-condensing  <b>Module Storage Conditions</b> Storage Temperature ..... -40-158°F (-40-70°C) Storage Humidity ..... 15-95%, non-condensing Exposure to Pressure ..... 500-1060 mbar

## Sensor Specifications

Application Site ..... Forehead	Ground Electrode ..... CB
Active Channels ..... 4	Reference Electrode ..... CT
Active Electrodes ..... L1, L2, R1, and R2	Duration of Use ..... Maximum of 24 hours
	Latex Content ..... Does not contain natural rubber latex

<sup>1</sup> Retrospective analysis of clinical data on file. <sup>2</sup> Purdon P et al. *Brit J of Anaesth*. 10.1093 46-57. <sup>3</sup> Narasway et al. *Critical Care Med*. 2002 Jul;30(7):1483-7. <sup>4</sup> Lobo, Francisco A., and Stefan Schraag. Limitations of anaesthesia depth monitoring. *Current Opinion in Anesthesiology*. 24, no. 6 (2011): 657-664. <sup>5</sup> Masimo Data on File.

The RD SedLine EEG Sensor is not licensed for sale in Canada.

For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, and precautions.

**Masimo U.S.**  
Tel: 1 877 4 Masimo  
info-america@masimo.com

**Masimo International**  
Tel: +41 32 720 1111  
info-international@masimo.com

